Prueba **RT-PCR**

Prepacol
S.A.S

¿Qué es RT- PCR?

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GeneFinder™ COVID-19 PLUS RealAmp kit used in association with ELITE InGenius® – Code: IFMR-45

A: Intended use

The GeneFinder™ COVID-19 PLUS RealAmp Kit is used for the Qualitative detection of SARS-CoV-2 virus, COVID-19 disease, through Reverse Transcription and Real-Time Polymerase Chain Reaction from RNA. The assay is compatible in combination with the instrument ELITE InGenius®


B: Amplified sequence

	Gene	Fluorophore
Target #1	RdRp gene	FAM
Target #2	E gene	Texas Red
Target #3	N gene	JOE
Internal Control (endogenous)	Rnase P gene	Cy5

C. Matrixes

› Alveolar Lavage Fluid (BAL)	Nasopharyngeal swabs (NPS)	Sputum
-------------------------------	----------------------------	--------

D. Kit content

COVID-19 PLUS Reaction Mixture (Purple)	COVID-19 PLUS Probe Mixture (Brown)	COVID-19 PLUS Positive Control (Red)	COVID-19 PLUS Negative Control (Green)	Maximum Shelf-life: 12 Months
				<p>Storage Temperature: -20°C</p>
<p>Reaction Mixture 1 tube of 1,050 µL 85 reactions per kit 10 freeze-thaw cycles</p>	<p>Probe Mixture 1 tube of 550 µL 85 reactions per kit 10 freeze-thaw cycles</p>	<p>Positive Control 1 tube of 50 µL 10 freeze-thaw cycles 8 tests per Kit</p>	<p>Negative Control 1 tube of 50 µL 10 freeze-thaw cycles 8 tests per 1 Kit</p>	<p>Transfer the controls in a 0,5 mL Sarstedt tube (Ref. 72730005)</p>

E. Master Mix Preparation

Reagent	1 Test	Samples N.	Total reactions
COVID-19 PLUS Reaction Mixture	10 µL	1 to 4	N + 1 rxn
COVID-19 PLUS Probe Mixture	5 µL	5 to 8	N + 2 rxn
Total volume	15 µL	9 to 12	N + 3 rxn

**Prepare the Mastermix
in a 2 mL Sarstedt tube
(Ref. 72694005)**

F. Material required not provided in the kit

ELITE InGenius instrument:	INT030
ELITE InGenius SP200 Extraction Cartridge:	INT032SP200
ELITE InGenius PCR Cassette:	INT035PCR
ELITE InGenius SP200 Consumable Set:	INT032CS
ELITE InGenius Waste Box:	F2102-000
Filter Tips 300:	TF-350-L-R-S
Sarstedt 0,5 mL Microtube, PP	72730005
Sarstedt 2 mL Microtube, PP	72694005

G. ELITe InGenius protocol

Sample volume	200 μ L
Total eluate volume	100 μ L
PCR eluate input volume	5 μ L
Complete PCR Mix volume	15 μ L
Frequency of controls	15 days

H. Result Interpretation

RdRP	E	N	Result
+	-	+	SARS-CoV-2 positive
+	+	-	
+	+	+	
+	-	-	Repeat to confirm Positive
	+	+	
	-	+	
	+		SARBeCoV

I. Performance

LoD (PCR)

10 copies (RdRp, E, N gene)

Reproducibility

CV < 5%

Repeatability

CV < 5%

J. Tube type collection

Copan Ref	Description
608CS01R	eNAT Swab*
359C	UTM Swab

*The eNAT medium contains guanidine thiocyanate

Note: Not use the Swab tube as primary tube but transfer 200 µL of sample into sonicator tube.

K. Sample Pretreatment (optional)

Sample inactivation is not required as shown in WHO guideline. The sample can be pretreated with guanidine thiocyanate or other lysis buffers under a Biosafety cabinet of class II (BSC2)

- Dispense 200 μ L of sample into the sonicator tubes
- Add 150 μ L of buffer containing denaturants
- Cover the sonicator tube with the cap
- Load the sample into the instrument removing the caps

Note: This procedure needs to be validated as per WHO guidelines

L. Procedures

The user is guided step-by-step by the ELITE InGenius software to prepare the run. All the steps: extraction, amplification and result interpretation are automatically performed. Three operational modes are available: complete run, or extraction only, or PCR only.

Before analysis

1. Switch on ELITE InGenius Identification with username and password Select the mode "Open"	Verify controls: COVID positive and negative controls in the "Control menu" NB: Both had to be run, approved and not expired	2. Thaw the COVID-19 PLUS reagent tubes. Vortex gently and spin down 5 secs Prepare the complete Master Mix as shown in the section E
--	--	--

Procedure 1 – Complete run: Extraction + PCR

1. Select “Perform Run” on the touch screen	2. Verify the extraction volumes: Input: “200 µL”, eluate: “100 µL”	3. Scan the sample barcodes with handheld barcode reader or type the sample ID
4. Select the “Assay protocol” of interest	5. Select the sample position: Primary tube or sonication tube	6. Load the complete Master Mix in the inventory block
7. Load: PCR cassette, Extraction cartridge, Elution tube, Tip and sample collected into sonicator tubes. For Swab use see J For inactivation step see K	8. Close the door Start the run	9. View, approve and store the results For result interpretation see H

Procedure 2 – PCR only

1 to 4 : Follow the Complete Run procedure described above	5. Select the protocol “PCR only” and set the sample position “Extra tube”	6. Load the extracted nucleic acid tubes in the rack n°4
7. Load the PCR cassette rack Load the Q-PCR Mix in the inventory block	8. Close the door Start the run	9. View, approve and store the results

Procedure 3 – Extraction only

1 to 4 : Follow the Complete Run procedure described above	5. Select the protocol “Extraction Only” and set the sample position: Secondary tube	6. Load: Extraction cartridge, Elution tube, Tip cassette, sonication tube and primary sample racks
8. Close the door Start the run	9. Archive the eluate sample	

Prepacol
S.A.S